


# Section 28: Bone Plate Analysis and Design

# Biomechanics of Bone Fx

- Weakest in Tension, Strongest in Compression
- Pure **Bending** → **Transverse Fx**
- **Torsion** → **Spiral Fx**
- **Shear** → **Oblique Fx**
- **Butterfly** due to **Bend + Shear**




# Screw fixation



- Rotary forces → compression between objects (inclined plane on spiral pulls object toward head)
- Four part construction: ***head, shaft, thread, tip***
- Thread defined by ***root diameter, thread diameter, pitch***

# Screws

- **Larger core diameter** has **higher resistance to fatigue & shear failure**
  - 4th power of the diameter
- **Pullout strength** (maximum force screw can support along its axis)
  - outer diameter, length of engagement, shear strength/density of bone






# Screw fixation




# Interfragmentary compression screw

- Lag screw can convert ***torque forces*** to a ***compressive force***
- Screw should be ***perpendicular*** to fracture preventing sliding of fragments when compressed
- ***Gliding hole*** => cortex under screw head drilled to thread diameter
- ***Thread hole*** => opposite glide hole; drilled to the core diameter of screw



# Interfragmentary compression screw

- This alone is well suited for avulsion, epiphyseal, metaphyseal and intraarticular fractures- absolute stability, no callus






# Plate and screw fixation

- For transverse or short oblique fx, screws must be combined with other internal fixation
- Principle based on converting *tension force* to ***compression force***
- Accomplished by placing plate on ***tension*** or ***convex*** side of bone

# Plate and screw fixation



28-10

From: Justice

# Plate and screw fixation

- Plates offer benefits of ***anatomic reduction*** and stability for early motion, but must be ***protected from early weight bearing***.
- Ideally ***3-4 screws*** on each side of fracture
- ***Prebending*** the plate to increase compression


# Biomechanics of Plate Fixation

- Function of the plate
  - Internal splint
  - Compression
- “The bone protects the plate”






# Biomechanics of Plate Fixation

- Fracture Gap /Comminution
  - Allows bending of plate with applied loads
  - Fatigue failure




# Plates & Bending failure

- Leaving gap opposite plate makes it a fulcrum
- increased stress at holes
- avoid holes over fracture sites
- greater the span between screws
  - less stiff
  - more bending




# Biomechanics of Plate Fixation

- Bone-Screw-Plate Relationship
  - Bone via compression
  - Plate via bone-plate friction
  - Screw via resistance to bending and pull out.



# Biomechanics of Plate Fixation


- Plates:
  - Bending stiffness proportional to the thickness (h) of the plate to the 3<sup>rd</sup> power.





# Biomechanics of Plate Fixation

- The screws closest to the fracture see the most forces.
- The construct rigidity decreases as the distance between the innermost screws increases.





**Screw Axial Force**

# Biomechanics of Plate Fixation

- Number of screws (cortices) recommended on each side of the fracture:

Forearm	3	(5-6)
Humerus	3-4	(6-8)
Tibia	4	(7-8)
Femur	4-5	(8)



# Prebending




28-19

From: Justice

- Dynamic Compression



- Tension Band Plate



**Fig. 3.2.2-23:** Tension band principle at the femur. The intact femur (a) is an eccentrically loaded bone with distraction or tensile forces laterally and compression on the medial side. In case of a fracture (b) the lateral fracture gap will open, whereas the medial will be compressed. If a plate is placed alongside the linea aspera of the femur (c), it will be under tension when loaded, thereby compressing the fracture gap, provided there is bony contact medially. If the plate is placed to the compression side (d), it is not able to prevent opening of the lateral gap (instability). If the medial cortex is not intact (e), the tension band principle cannot work, because of lack of buttress (see [chapter 3.2.3](#)).